

IATROSCAN INSTRUMENT APPLICATION

No. 16

Analysis of Triglyceride Molecular Species using the New SIII  
CHROMARODS (silver nitrate impregnated)

Analysis of Triglyceride Molecular Species  
using Silver Nitrate impregnated Chromarods (new type SIII)

Triglyceride molecular species can easily be separated and analysed quantitatively on the basis of different degrees of unsaturation of the constituent fatty acids by the IATROSCAN TLC/FID using the new type CHROMAROD SIII treated with silver nitrate. This analysis can be used as a simple method of detection for oils and fats.

1. Chromarods treated with silver nitrate.

A. Preparation of the CHROMARODS:

Make up a 3% silver nitrate Acetonitrile solution:


The silver nitrate (0.9g), fully dissolved in 30 ml of acetonitrile, is placed in a test tube. Cover the test tube with aluminium foil to exclude the light.

The preparation of a silver nitrate-treated CHROMAROD is then carried out in the following manner:

- ① CHROMAROD-SIII is first subjected to a blank-scan in the MK IV IATROSCAN to remove any absorbed organic substances.
- ② Immerse the CHROMAROD into 3% silver nitrate Acetonitrile solution for 5 minutes.
- ③ The CHROMAROD is then removed from the solution and dried in an oven at 120 degrees C for 10 minutes.

In this way the CHROMAROD is properly treated.

Fig 1 The preparation sequence for CHROMARODS treated with silver nitrate



B. Washing and regeneration of the CHROMARODS treated with silver nitrate:

Once CHROMARODS treated with silver nitrate are used, the effectiveness of the silver nitrate is reduced. The following procedure for washing and regeneration is, therefore, recommended.

- ① Immerse CHROMARODS in concentrated nitric acid overnight.
- ② Prepare two test tubes filled with pure water.
- ③ Remove CHROMARODS from the nitric acid and place them in the test tubes ② in order to wash them.
- ④ The CHROMARODS are then dried in an oven at 120 degrees C for one hour.
- ⑤ Treat CHROMARODS with silver nitrate in same manner as stated in A, ① through ③

The RODS can be repeatedly used after being washed and reactivated according to the above procedures (repeated use for up to ten times or more does not adversely affect the separation performance).

2. Analytical Method

- ① Spot a chloroform solution (1 µl) containing 10 to 20 µg of lipid onto the silver nitrate-treated CHROMARODS prepared according to Section 1.

- ② Develop the spot using

Benzene : Ethylether : Formic acid  
68            2            0.1            ..... (solvent mixture A)

Benzene : Ethylether : Formic acid  
65            5            0.1            ..... (solvent mixture B)

Benzene : Chloroform : Acetic acid  
63            7            0.7            ..... (solvent mixture C)

The solvent mixture A is suitable for the analysis of samples containing constituents having up to about 5 double combinations with a lesser degree of unsaturation. Use solvent mixtures B and C for relatively higher degrees of unsaturation.


Remarks: While developing the spot, the development chamber is required to be wrapped in aluminium foil to prevent decomposition of silver nitrate by exposure to light.

- ③ Dry the CHROMARODS after development at 120 degrees C in a dryer for 5 minutes to remove all solvent.
- ④ Measure the CHROMARODS in the IATROSCAN analyser after solvent removal.

Conforming to the analytical procedures explained above, a standard mixture of triglyceride (tripalmitin, triolein and trilinolein) developed in solvent mixtures A, B and C exhibit chromatographic separations shown in Fig.2. Identification of oils and fats is also shown in Figs. 3 and 4 respectively.


Fig. 2 Chromatograms obtained with a mixture of standard Triglycerides by using silver nitrate-treated CHROMARODS

Fig 2-1. Solvent mixture A




CAL. METHOD 00					
SF .100000e+03 PA .100000e+01 PB .100000e+01					
NO.	NAME	RT	A OR H	MK	CONC
1	P3	0.173	12064		36.2988
2	O3	0.288	12254		36.8621
3	L3	0.461	8918	H	26.8320
TOTAL			33236		100.0000

Fig 2-2. Solvent mixture B



CAL. METHOD 00					
SF .100000e+03 PA .100000e+01 PB .100000e+01					
NO.	NAME	RT	A OR H	MK	CONC
1	P3	0.162	26339	H	34.6696
2	O3	0.204	26868		35.3670
3		0.276	899	H	1.1311
4	L3	0.459	21913		28.8402
TOTAL			75922		100.0000

Fig 2-3. Solvent mixture C



CAL. METHOD 00					
SF .100000e+03 PA .100000e+01 PB .100000e+01					
NO.	NAME	RT	A OR H	MK	CONC
1	P3	0.187	24530		35.7330
2	O3	0.258	23836		34.7222
3		0.328	148		0.2166
4		0.384	627		0.3144
5	L3	0.457	19438	H	28.3153
6		0.492	66	H	0.0973
TOTAL			66648		100.0000

Identifications :

P3 : Tripalmitin

O3 : Triolein

L3 : Trilinolein

Mobile phase :

Solvent mixture A

Benzene : Ethyl ether : Formic acid

68 : 2 : 0.1

Solvent mixture B

Benzene : Ethyl ether : Formic acid

65 : 5 : 0.1

Solvent mixture C

Benzene : Chloroform : Acetic acid

63 : 7 : 0.7

CONDITIONS :

Stationary phase : CHROMAROD-SIII (3% silver nitrate impregnated)

Gas flow : H<sub>2</sub> 160ml/min, Air 2.0l/min

Scanning speed : 30sec/scan

Introcorder TC-11

Playback attenuation : 16-32mv f.s


3. Application for identification of oils and fats. Figs. 3 and 4 show chromatograms obtained with samples of fats, palmoil, cocoa butter, coconut oil and beef tallow.

Fig. 3 Patterns of separated components of fats using CHROMARODS treated with silver nitrate.

Fig 3-1. Palm oil

Solvent mixture A

SAMPLE 8 15156 JUNE 30 1987




CAL. METHOD 00						
SF PA PF						
.100000e+03 .100000e+01 .100000e+01						
NO.	NAME	RT	A OR H	NK	CONC	
1		0.216	17763	H	35.9472	
2		0.261	13152		27.2291	
3		0.311	8085	H	16.7402	
4		0.384	5425	"	11.2302	
5		0.446	4224	"	8.7459	
TOTAL			48301		100.0000	

Fig 3-2. Cacao oil

Solvent mixture A

SAMPLE 9 15159 JUNE 30 1987




CAL. METHOD 00						
SF PA PF						
.100000e+03 .100000e+01 .100000e+01						
NO.	NAME	RT	A OR H	NK	CONC	
1		0.210	46371	H	78.5560	
2		0.248	8922		15.1149	
3		0.324	1794	H	3.0223	
4		0.416	1751		3.3067	
TOTAL			59029		100.0000	

Fig 3-3. Coconut oil

Solvent mixture A

PLAY BACK ATTENUATE 16  
SAMPLE 10 16101 JUNE 30 1987




CAL. METHOD 00						
SF PA PF						
.100000e+03 .100000e+01 .100000e+01						
NO.	NAME	RT	A OR H	NK	CONC	
1		0.205	31367	H	30.8494	
2		0.274	1242		3.5280	
3		0.332	576	H	1.6710	
4		0.363	625	H	2.7207	
5		0.396	514	H	1.4911	
TOTAL			34526		100.0000	

Fig 3-4. Beef tallow

Solvent mixture A

SAMPLE 23 16133 JUNE 30 1987




CAL. METHOD 00						
SF PA PF						
.100000e+03 .100000e+01 .100000e+01						
NO.	NAME	RT	A OR H	NK	CONC	
1		0.155	448		.553	
2		0.213	2880		91.9426	
3		0.261	10993	H	13.7678	
4		0.331	2420	H	4.1327	
5		0.370	5891	H	8.4324	
6		0.427	3693	H	6.3056	
7		0.468	6029	H	10.2934	
TOTAL			56575		100.0000	

Fig. 4 Patterns of separated components of fats using CHROMARODS treated with silver nitrate.

Fig 4-1. Olive oil


Solvent mixture A



NO.	NAME	RT	A OR H	HK	CONC
1		0.207	2661	H	3.8898
2		0.247	18369	H	26.8474
3		0.296	34707	H	50.7258
4		0.378	8455	H	12.3580
5		0.451	4227	H	5.1788
TOTAL			68421		100.0000

Fig 4-2. Soybean oil


Solvent mixture B



NO.	NAME	RT	A OR H	HK	CONC
1		0.213	17666	H	28.3136
2		0.238	18640	H	21.4340
3		0.271	14445	H	16.6108
4		0.310	30207	H	23.2355
5		0.382	4665	H	5.3645
6		0.411	8537	H	9.8168
7		0.465	2804		3.2246
TOTAL			86267		100.0000

Fig 4-3. Sunflower oil


Solvent mixture B



NO.	NAME	RT	A OR H	HK	CONC
1		0.196	2275	H	3.5456
2		0.211	3782	H	9.3844
3		0.236	10906	H	24.8248
4		0.276	10192	H	25.2627
5		0.324	13412	H	33.2764
6		0.449	687		1.7058
TOTAL			40306		100.0000

Fig 4-4. Safflower oil

Solvent mixture C




NO.	NAME	RT	A OR H	HK	CONC
1		0.200	220	H	0.4282
2		0.224	1432	H	2.7805
3		0.258	3267	H	6.3430
4		0.288	11732	H	21.9969
5		0.333	8182	H	13.8324
6		0.361	25505	H	49.5060
7		0.457	1577	H	3.0627
TOTAL			51519		100.0000

Fig 4-5. Linseed oil

Solvent mixture B

SAMPLE 14 14:20 JULY 01 1987




NO.	NAME	RT	A OR H	HK	CONC
1		0.121	1405	H	3.0537
2		0.216	1306	H	2.8385
3		0.253	1453	H	3.1581
4		0.287	2931	H	6.2625
5		0.326	5254	H	11.4414
6		0.358	6354	H	18.1799
7		0.462	25336	H	55.0458
TOTAL			46910		100.0000

Fig 4-6. Tung oil

Solvent mixture B

SAMPLE 18 14:25 JULY 01 1987




NO.	NAME	RT	A OR H	HK	CONC
1		0.203	13548	H	38.5773
2		0.241	5436	H	15.4810
3		0.306	5900	H	11.1072
4		0.369	5222	H	14.8712
5		0.460	7911	H	19.9531
TOTAL			35119		100.0000

Fig 4-7. Rape seed oil

Solvent mixture B

SAMPLE 6 15:52 JUNE 30 1987



NO.	NAME	RT	A OR H	HK	CONC
1		0.172	1521	H	2.9324
2		0.258	3876	H	7.2880
3		0.308	15637	H	29.5119
4		0.362	11492	H	21.5984
5		0.434	20541	H	38.6191
TOTAL			53109		100.0000

CONDITIONS :

Stationary phase

CHROMAROD-SIII (3% AgNO<sub>3</sub> impregnated)

Mobile phase :

Solvent mixture A

Benzene : Ethyl ether : Formic acid

68 : 2 : 0.1

Solvent mixture B

Benzene : Ethyl ether : Formic acid

65 : 5 : 0.1

Solvent mixture C

Benzene : Chloroform : Acetic acid

63 : 7 : 0.7

Scanning speed : 30sec/min

Iatrocoder TC-11

Playback attenuation : 16

Reference :

- 1) M. Tanaka, T. Itoh and H. Kaneko : Yukagaku 28(2):96, 1979.
- 2) Sebedio, J.L., Ackman, R.G. : J. Chromatogr. Sci. 19 : 552, 1981.

As shown in Figs. 3 and 4, solid lipids contain numbers of triglyceride molecular species of a relatively lesser degree of unsaturation. On the other hand, liquid lipids contain numbers of triglyceride molecular species of a higher degree of unsaturation which reflect the composition ratio of the constituents of fatty acids respectively.

Depending on the type of oils and fats, specific chromatograms of triglyceride constituents can be obtained and are applicable for the identification of oils and fats.

#### 4. Analytical results



Table 1 shows reproducibility relating to a standard triglyceride analysed using a CHROMAROD-SIII treated with 3% silver nitrate, and the area percentage obtained. When the total amount spotted varies from 4 to 20  $\mu$ g in the case of olive oil, the relationship between total weight and response can be seen (Fig. 5).

Table 1 - Reproducibility of percentage of standard triglycerides.

Rod No.	Tripalmitin %	Triolein %	Trilinolein %
1	34.6	37.6	27.8
2	34.9	36.8	28.3
3	33.7	36.8	29.5
4	34.9	36.4	28.7
5	34.9	36.9	28.2
6	34.6	36.7	28.7
7	34.9	37.1	28.0
8	35.4	37.5	27.1
9	36.3	36.9	26.8
10	35.0	36.7	28.3
$\bar{X}$	34.9	36.9	28.2
CV%	1.9	0.8	2.8



Fig. 5 - Relationships between the amount of olive oil spotted and the response



Stationary phase:

CHROMAROD-SIII (3% AgNO<sub>3</sub> impregnated)

Mobile phase:

Solvent mixture A

Benzene: Ethyl ether : Formica acid  
68            2                    0.1